



**DIAG**  
IMPORT & EXPORT

Prices exclude VAT & Discount

**DIAG MAY/JUNE SPECIALS!**

Tel: 011 794 1577

Fax: 011 794 1852

Email: orders@diagsa.co.za

**Products displayed less 10% - While Stocks Last or 30 June 2017.**

Valid from 15 May 2017 – 30 June 2017

**SKIPPER CARRIERS:**

MAR0001	1F	R231.27 excl
MAR0002	2F	R328.37 excl
MAR0004	3F	R420.89 excl



**FURBED'S:**

LO91581	150X150	R231.87 excl
LO91518	100X150	R161.42 excl
LO91579	100X100	R113.23 excl
LO91584	75X100	R80.55 excl
LO91582	65X75	R52.77 excl
LO91580	50X65	R33.76 excl

**OSTER GOLDEN A5 CLIPPER:**

Oster Golden A5 Clipper	LO91585	R3389.23 excl
-------------------------	---------	---------------



**B BRAUN SUTURE MATERIAL: [Please email for Prices](#)**

**WIDE RANGE AVAILABLE IN DIFFERENT SIZES AND NEEDLES.**

**Monosyn**

Monosyn is a mid-term absorbable monofilament and uncoated suture material made from glyconate. Retention of 50% knot tensile strength after 14 days, complete absorption after 60 – 90 days.



**Monoplus**

Monoplus is a long-term absorbable monofilament and uncoated suture made of polydioxanone. Retention of 50% knot tensile strength in 35 days.



**Softcat**

Softcat Chrom is a natural, twisted absorbable catgut suture made of pure collagen and treated with chromic salts. Retention of 50% knot tensile strength in 14 – 21 days. Complete absorption in 90 days.



**Dafilon**

Dafilon is a non-absorbable nylon monofilament suture.



**Novosyn**

Novosyn is a mid-term absorbable braided and coated synthetic suture made of Poly (glycolide-co-l-lactid 90/10) coated with Poly (glycolide-co-l-lactide 35/65) + CaSt Retention of 50% knot tensile strength after 21 days. Complete mass absorption: 56 up to 70 days.



**Suture Reels**

Chromic Catgut  
Dafilon Blue  
Safil Violet

